


SCRATCHED WITH A YELLOW TONE

Visible material, more assured




THREE CLEANING MODES TO CAREFOR EVERY NEED

Different dental cleaning modes, master free



Soft

For beginners, gently remove tartar and stains



Comfort

Removes stubborn tartar such as tobacco stains, coffee stains,tea stains,etc.




Strong

Remove formed calculus and clean teeth more powerfully



HOW MUCH DO YOU KNOW THE DANGERS OF DENTAL PROBLEMS

One machine to deal with all kinds of dental problems




Red, swollen and painful gums, always bleeding from brushing



Stubborn smoke stains, tea

stains, yellow and black teeth



Tartar calculus,accompanied by bad breath



Plaque accumulates, causing tooth decay

WE REDEFINE, VISUAL ULTRSONIC DENTAL SCALER

High-efficiency ultrasonic, high-power resonant tooth cleaning, with 2 million high-definition cameras, rea izes visual footh cleaning, and oral stains are clearly visible, and tooth dirt can be accurately cleaned.




MAGNETIC LEVITATION SONIC MOTOR, REFRESH AND ENJOY CLEANNESS BETWEEN TEETH

The magnetic suspension high-speed motor produces high-frequency sound waves that oscillate in a wide range, forming a cleaning force with poothpaste foam for deep cleaning.

VISIBLE FOUR-AXIS GYROSCOPE 360°NO DEAD ANGLE PRECISE TEETH CLEANING

Built-in four-axis directional gyroscope, the camera's angle of view remains unchanged when the body rotates, and 360° accurate tooth cleaning, giving you a safe and comfortable tooth cleaning experience.



2 MILLION PROFESSIONAL CAMERAS, VISUALLY CLEANING TEETH CLEANER

The camera is equipped with macro processing technology and built-in multiple sensors, allowing you to enjoy the 2 million high-definition pixel experience





before use


Teeth turn yellow and black, stubborn tartar is obvious

After use

Dental plaque and tartar are significantly reduced


SLIM CLEANING HEAD FOR PRECISE CLEANSING DEEP INTO THE TOOTH FRONT


The high-definition mobile phone APP shows that itgoes deep into the teeth to fully clean the dent alcalculus. There is nowhere to hide the dental plaque. It turns into a deep cleaning between the teeth.



REFUSAL TO BE AFRAID, INTELLIGENT RECOGNITION DOES NOT HURT THE GUMS

It will only work when it touches teeth or hard objects, and it will automatically stop working when it touches gums or soft objects.





Touching teeth

Touching gums

Ultrasonic vibration, start to work

Stop working and protect your gums

SMALL AND PORTABLE LIFE WATERPROOF

The head of the fuselage has a waterproof design and is highly waterproof. After cleaning the tartar, the lens can be washed to prevent the growth of bacteria.



Parameter:

Product name: Visual ultrasonic dental scaler

Pixel:2 million

Charging time:about 2-3 hours

Use time:about 120 minutes of continuous use