



Tangshan UMG Medical Instrument Co., Ltd

ADD: Extraordinary Base, No. 103 Jianshe North Road, High-tech Zone, Tangshan City,
Hebei Province
TEL:+86-315-6314398 +86-186-33395498 FAX:+86-315-6314577
EMAIL:info@umgroupen.com <https://www.umgroupen.com>

Dental Rotary Instrument



Professional manufacturer
of Dental Tool
Edition 2025

High Gloss Polishing RA Kit

Multi-function Rubber Wheel Kit

Order No: HP 3403G

Packing: 3 pcs/kit

For High Brightness Polishing of Composite Resin



Introduction

Characteristic : special for composite material and resin polishing , fast polishing speed , low heat generation , wonderful high light effect

- Steps:
- 1) Blue Polisher : Eliminating surfacescratches
 - 2) Brown Polisher: Pre-Polishing
 - 3) Grey Polisher : High-gloss polishing
- (Recommend R.P.M7000-12000/Min Max.15000/Min)

Caution : Do not exceed the recommended maximum speed, avoid sharp corners and slit polishing, and prohibit reverse use.

Composite Rubber Wheel Kit

Order No: RA 3306G

Packing: 6 pcs/kit

For High Brightness Polishing of Composite Resin



Introduction

Characteristic : special for High Brightness Polishing of Composite Resin, Zirconia, Porcelain, Metal, fast polishing speed , low heat generation , wonderful high light effect

- Steps:
- 1) Blue Polisher : Eliminating surfacescratches
 - 2) Brown Polisher: Pre-Polishing
 - 3) Grey Polisher : High-gloss polishing
- (Recommend R.P.M7000-12000/Min Max.15000/Min)

Caution : Do not exceed the recommended maximum speed, avoid sharp corners and slit polishing, and prohibit reverse use.

Multi-function Rubber Wheel Kit

Order No: RA 3403G

Packing: 3 pcs/kit

For High Brightness Polishing of Composite Resin,Zirconia,Porcelain,Metal



Introduction

Characteristic: Special for composite material and resin polishing, fast polishing speed, low heat generation, wonderful high light effect

- Steps:
- 1) Blue polisher: Eliminating surface scratches
 - 2) Brown polisher: Pre-polishing
 - 3) Grey Polisher: High gloss polishing
- (Recommend R.P.M7000-12000/Min Max. 15000/min)


High Gloss Polishing RA Kit

Polishing Kit For Glass Ceramic

Order No: RA 1212 D

Packing : 12 pcs/kit

Item No.: CD 3144, CD 3064, RD 3114, RD 3144, RD 3414, RD 3204, RD 3117A, RD 3147A , RD 3417A, RD 3207A



Introduction

Characteristic: Fast grinding, fast heat dissipation, excellent polishing effect.


- Steps:
- 1) Ceramic grinders are used for grinding the Maxillofacial.
 - 2) Grey polishers are used for smoothing and pre-polishing.
 - 3) Light pink polishers are used for high gloss polishing.
- Effect:
- 1) The sharpness of ceramic diamond grinder can leveling the connection channel between crown and porcelain quickly.
 - 2) Low heat production can avoid the damage to restoration because of overheating.
 - 3) Shorten the operating time beside the chair and immediately achievingthe same high gloss polishing effect as after the glaze.
- Caution: Control the pressure and rotate speed when grinding the connection channel between teeth and porcelain toavoid the damage of restoration. Pre-polishing is the base of high gloss polishing. So Pre-polishing is indispensable.

Polishing Kit For Enamic

Order No: RA 1312 D

Packing : 12 pcs/kit

Item No.: RD 3113A, RD 3143A, RD 3203A, RD 3413A, RD 3114A, RD 3144A, RD 3414A, RD 3204A, RD 3117B, RD 3147B, RD 3417B, RD 3207B



Introduction

Characteristic: Fast grinding, long service life, low heat,production, good polishing effect.

- Steps:
- 1) Grey polishers are used for grinding the Maxillofacial.
 - 2) Deep pink polishers are used for pre-polishing.
 - 3) Greyish- white polishers are used for high gloss polishing.
- Effect:
- 1) Low heat production can avoid the damage to restoration because of overheating.
 - 2) By adjusting the rotate speed and pressure of deeppink polisher, achieving smoothing and polishing effect.
 - 3) Shorten the operating time beside the chair and immediately achievingthe same high gloss polishing effect as after the glaze.
- Caution: Control the pressure and rotate speed when grinding the connection channel between teeth and porcelain to avoid the damage of restoration. Pre-polishing is the base of high gloss polishing. So Pre-polishing is indispensable.


High Gloss Polishing RA Kit

Polishing Kit For Zirconia

Order No: RA 0112 D

Packing(pcs): 12

Item No.: CD 3064 x2, CD 3144 x2, RD 3114, RD 3144, RD 3414, RD 3204, RD 3115, RD 3145 RD 3415, RD 3205



Introduction

Characteristic: Fast cutting, long service life ,low heat production,good polishing effect.

(No water is required!)

Steps: 1) Ceramic diamond grinders are used for grinding the Maxillofacial morphology.

2) Grey polishers are used for prepolishing.

3) Yellow polishers are used for high gloss polishing.

Effect: 1) Take advantage of the special properties of ceramic diamond grinders(without water conditions) to fast finish the crystal zirconia restorations, reduced the complicated operation process of the traditional diamond bur with the water conditionin.

2) Take advantage of the special properties of the rubber diamond to achieve the high gloss polishing effect.


Caution: Without too much pressure. According to the recommended speed,with better effect.

Polishing Kit For Porcelain

Order No: RA 0212 D

Packing : 12 pcs/kit

Item No.: CW 3064m x2, CW 3114m x2, SH 3113, SH 3143m,SH 3413m, SH 3203m,RD 3116,RD 3146,RD 3416,RD 3206



Introduction

Characteristic: Fast speed grinding, long service life, fast heat dissipation, excellent polishing effect.

Steps: 1) White stones are used for grinding the maxillofacial.

2) Pink silicone polishers are used for prepolishing.

3) White rubber diamond polishers are used for high gloss polishing.

Effect: 1) Exquisite pure alumina sand decreased the vibration to the porcelain teeth occlusion and adjacent surface finishing, avoid carbon pollution at the same time.

2) Good heat dissipation performance of silicon can avoid the porcelain cracked.

3) Polishers with diamond grains have efficient polishing effect, wear-resistant durable.

4) High gloss polishing effect can avoid secondary rework of the color deviation and porcelain brittle risk

Caution: Polishing the near edge of the gold crown porcelain junction, please note that the pressure of the hand, grinding direction, speed, avoid collapse porcelain, porcelain crack and crack phenomenon. Need to prepolishing,lay a foundation for high gloss polishing.


High Gloss Polishing RA Kit

Polishing Kit For Composite

Order No: RA 0312 D

Packing : 12 pcs/kit

Item No.: CW 3064m x2,CW 3114m x2,SH 3114m,SH 3144m,SH 3414m, SH 3204m,RD 3117A,RD 3147A,RD 3417A,RD 3207A



Introduction

Characteristic: Fast speed grinding, long service life, fast heat dissipation, excellent polishing effect.

Steps: 1) White stones are used for adjusting anatomical shape and remove excesscomposite to finish shaping.

2) Light blue polishers are used for eliminating surfacescratches and prepolishing.

3) Light pink polishers are used for high gloss polishingof the surface.

Effect: 1) Fast finishing the balance composite after light cure.

2) After polishing, without any micro leakage at the contact point and surface of composite with the tooth.

3) Experimental results show that the perfect polishing effect can prolong the service life of 1-2 years.

Caution: When finishing and polishing, control the pressure and rotate speed, avoid damage the enamel.

Polishing Kit For Composite

Order No: RA0309

Packing : 9pcs/kit

Type: CW3064, CW3114, CW3514, SH3133, SH3203, SH3413, SH3134, SH3204, SH3414




Usage: Special for porcelain and resin teeth. All silicone polishers for smoothing and initial polishing easily.

Interproximal enamel reduction kits

Order No: FG 3205 E

Packing : 5pcs/kit

Item No: G/TC-07F, G/TC-08F, G/TC-08EF, G/TR-51F, G/TR-51EF



High Gloss Polishing RA Kit

Universal polishing kit

Order No: RA 0406E
Packing: 6pcs/kit
Item No: CD 3064, RD 314 4B, RD 3145A, RD 3147A, RD3377C



RA0406E (high speed, low heatness, high gloss)

01. Steps for Zirconia polishing

- 1) CD3064 light yellow ceramic diamond polisher is used for preliminary grinding and finishing.
- 2) RD3144B blue polisher is used for preliminary polishing.
- 3) RD3145A yellow polisher for high-brightness polishing.

Pay attention: Controlling the pressure and speed during polishing, and pre-polishing is required to lay the foundation for bright polishing.

Working diameter (mm)	3x7	4.5x9	4.5x9
Grit	Fine	Fine	Extra Fine
Package(pcs)	5/100	5/100	5/100
Model	CD 3064	RD 314 4B	RD 3145A
step	First	Second	Last

02. Polishing steps for glass ceramics

- 1) CD3064 light yellow ceramic diamond polisher is used for preliminary grinding and finishing.
- 2) RD3144B blue polisher is used for preliminary polishing.
- 3) RD3147A light pink polisher is used for high-brightness polishing.
- 4) RD3377C gray white polisher is used for high brightness polishing.

Pay attention: It is not necessary to use too much pressure during polishing. Use the recommended rotating speed, and do a good job of pre-polishing to avoid damaging the mouth due to high temperature generated by high speed, so as to lay the foundation for bright polishing.

Working diameter (mm)	3x7	4.5x9	4.5x9	11x1.6
Grit	Fine	Fine	Extra Fine	Extra Fine
Package(pcs)	5/100	5/100	5/100	5/100
Model	CD 3064	RD 314 4B	RD 3147A	RD337 7C
step	First	Second	Last	Available

03. Steps for resin polishing

- 1) CD3064 light yellow ceramic diamond polisher is used for preliminary grinding and finishing.
- 2) RD3144B blue polisher is used for preliminary polishing.
- 3) RD3147A light pink polisher is used for high-brightness polishing.
- 4) RD3377C gray white polisher is used for high brightness polishing.

Pay attention: During trimming and polishing, control the rotation speed to avoid damaging the enamel.

Working diameter (mm)	3x7	4.5x9	4.5x9	11x1.6
Grit	Fine	Fine	Extra Fine	Extra Fine
Package(pcs)	5/100	5/100	5/100	5/100
Model	CD 3064	RD 314 4B	RD 3147A	RD337 7C
step	First	Second	Last	Available

High Gloss Polishing HP Kit

Polishing Kit For Enamic

Order No: HP 1310 D
Packing : 10 pcs/kit
Item No.: TH 2013, RD 2133A, RD 2373A, RD 2423A, RD 2134A, RD 2374A, RD 2424A, RD 2137B, RD 2377B, RD 2427B



Introduction

Characteristic: Fast grinding, long service life, low heat production, good polishing effect.

Steps: 1) Carbide burs are used for leveling the connection channel between the teeth and porcelain. (Note: Elastic ceramics contain resin component, do not choose other grinding burs)


- 2) Grey polishers are used for leveling the surface of the crown.
- 3) Dark pink polishers are used for pre-polishing.
- 4) Greyish- white polishers are used for high gloss polishing.

Effect: 1) The sharpness of carbide burs can leveling the connection channel between the teeth and porcelain quickly.
2) Low heat production can avoid the damage of restoration because of overheating.
3) By adjusting the rotating speed and pressure of dark pink polisher, achieving leveling and polishing effect.
4) Shorten the operating time beside the chair and immediately achieving the same high gloss polishing effect as after the glaze.

Caution: Please Control the pressure and rotate speed when grinding the connection channel between teeth and porcelain to avoid the damage of restoration. Pre-polishing is the base of high gloss polishing. So Pre-polishing is indispensable .

Polishing Kit For Glass-Ceramics

Order No: HP 1210 D
Packing : 10 pcs/kit
Item No.: CD 2124 , CD 2024, RD 2064 , RD 2134 , RD 2374 , RD 2424 , RD 2067A , RD 2137A , RD 2377A, RD 2427A



Introduction

Characteristic: Fast grinding, long service life, low heat production, good polishing effect.

Steps: 1) Ceramic grinders are used for leveling the connection Channel between the teeth and porcelain.

- 2) Grey polishers are used for pre-polishing.
- 3) Light pink polishers are used for high gloss polishing.


Effect: 1) The sharpness of ceramic diamond grinder can leveling the connection channel between teeth and porcelain quickly.
2) Low heat production can avoid the damage to restoration because of overheating.
3) Shorten the operating time beside the chair and immediately achieving the same high gloss polishing effect as after the glaze.

Caution: Please Control the pressure and rotate speed when grinding the connection channel between teeth and porcelain to avoid the damage of restoration. Pre-polishing is the base of high gloss polishing. So Pre-polishing is indispensable.

High Gloss Polishing HP Kit

Polishing Kit For Porcelain

Order No: HP 0209 D
Packing : 9 pcs/kit
Item No.: CW 2054m,CW 2134m, CW 2114,SH 2123m ,SH 2323m,SH 2413, RD 2126, RD 2326,RD 2416



Introduction

- Characteristic: Fast speed grinding, long service life, fast heat dissipation, excellent polishing effect.
- Steps:
- 1) White stones are used for grinding the maxillofacial morphology.
 - 2) Pink silicone polishers are used for smoothing and prepolishing.
 - 3) White rubber diamond polishers are used for high gloss polishing.
- Effect:
- 1) Exquisite pure alumina sand decreased the vibration to the porcelain teeth occlusion and adjacent surface finishing, avoid carbon pollution at the same time.
 - 2) Good heat dissipation performance of silicon can avoid the porcelain cracked.
 - 3) Polishers with diamond grains have efficient polishing effect, wear-resistant durable.
 - 4) High gloss polishing effect can avoid secondary rework of the color deviation and porcelain brittle risk
- Caution:
- Polishing the near edge of the gold crown porcelain junction, please note that the pressure of the hand, grinding direction, speed, avoid collapse porcelain, porcelain crack and crack phenomenon. Need to prepolishing, lay a foundation for high gloss polishing.

Polishing Kit For Zirconia

Order No: HP 0109 D
Packing : 9 pcs/kit
Item No.: CD 2124,CD 2024, CD 2324,RD 2124,RD 2324,RD 2414,RD 2125, RD 2325,RD 2415



Introduction

- Characteristic: Fast cutting, long service life, radiating faster, good polishing effect. (No water is required!)
- Steps:
- 1) Ceramic diamond grinders are used for grinding the Maxillofacial morphology.
 - 2) Grey polishers are used for leveling and prepolishing.
 - 3) Yellow polishers are used for high gloss polishing.
- Effect:
- 1) Take advantage of the special properties of ceramic diamond grinders(without water conditions) to fast finish the crystal zirconia restorations,reduced the complicated operation process of the traditional diamond bur with the water conditionin.
 - 2) Take advantage of the special properties of the rubber diamond to achievethe high gloss polishing effect.
- Caution:
- Without too much pressure. According to the recommended speed,with better effect.

High Gloss Polishing HP Kit

Polishing Kit For Acrylic

Order No: HP 0308 D
Packing : 8 pcs/kit
Item No.:TI 2064,TI 2074,TI 2053,SK 2135m,SK 2145m,SK 2025m,RD 2137, RD 2127



Introduction

- Characteristic: Fast cutting, long service life, radiating faster,good polishing effect.
- Steps:
- 1) Carbide burs are used for grinding shape.
 - 2) Yellow polishers are used for leveling and prepolishing.
 - 3) Pink polishers are used for high gloss polishing.
- Effect:
- 1) The carbide burs can grinding acrylic resin base lip,clip, lingual and occlusal surfaces lacing margin quickly
 - 2) High quality polishers can high gloss polishing the acrylic.
 - 3) Perfect polishing effect can improve the comfort for patient.
- Caution:
- Without too much pressure. According to the recommended speed,with better effect.

HP Polisher

HP Polisher

Order No:UM-HP001

Packing:25 models,
10pcs/board100pcs/box



RA Polisher

RA Polisher

Order No:UM-RA001


Packing:100pcs/bottle



HP Zirconia Ceramics Polisher

Order No:UM-HP002

Packing: 5 models, 1pc/box



RA Polisher for polishing light
curing and grinding teeth

Order No:UM-RA002

Packing:12pcs/set



HP Diamond Spiral Disc

Order No:UM-HP003

Packing:L,M,S 1pc/box



RA Diamond Spiral Disc

Order No:UM-RA003

Packing:silver handle 6pcs/box



RA Diamond Spiral Disc

Order No:UM-RA004
Packing:gold handle
6pcs/box



RA Spiral & eagle claw Wheel kit

Order No:UM-RA005
Packing:gold handle N2pcs spiral+2pcs eagle claw wheel/box



RA Resin polished eagle claw wheel kit

Order No:UM-RA006
Packing:gold handle, white, pink, blue 4pcs/box



RA Spiral & eagle claw Wheel kit

Order No:UM-RA007
Packing:gold handle 3pcs spiral+3pcs eagle claw wheel/box




Kit of Preparation

Polishing Kit For Composite

Order No: FG 0312 D

Packing : 12 pcs/kit

Item No.: BR-13,BR-13F,BR-15F,SR-33C,SR-33F,RA G1-014,FO-15F,TR-31F,FO-28F,MI-01,SH 3143m,RD 3115



Introduction


1. BR-13: Used for small and medium-sized hole corrosion removal
2. BR-13F: Used for the preparation of round concave edge of V-shaped hole
3. BR-15F: Used for preparation of labial concave shoulder edge of front teeth
4. SR-33C: For corrosion removal and hole edge trimming
5. SR-33F: For hole edge trimming and jaw edge adjustment
6. RA G1-014: For subsequent decay removal with texture and thoroughness
7. FO-15F: For surface texture processing
8. TR-31F: Used for trimming the adjacent surface and flash of Class II hole
9. FO-28F: Used for lingual and posterior jaw edge adjustment
10. MI-01: For removal of minor caries
11. SH 3143m: For polishing the cavity edge of natural teeth
12. RD 3115: For preliminary polishing after resin repair

Kit Of Preparation Anterior and Posterior Teeth for ceramics/ zirconia

Order No: FG 1112 D

Packing : 12 pcs/kit

Item No.: SO-20,BR-15,(Patent number:2014201580245),EX-54,EX-55,TR-33,TR-D10,FO-23,SF-D41F,TR-08F,TR-D36F,TR-D36F,SH 3414



Introduction

1. SO-20 : Open the collar near the distal, reduction is 1.0 mm, please do not hurt the teeth. (Use matrix band to protect the next teeth.)
2. BR-15 : Prepare tongue surface edge and position the shoulder 0.8mm,Depth is shallow concave.
3. EX-52 : For anterior biting edge positioning deep 1.5-2.0mm, it is easy to control the depth.For Posterior biting edge position deep 1.5-2.0mm, it is easy to control the depth.
4. EX-54 : For anterior labial surface transverse positioning deep 1.2-1.5mm, it is easy to control the depth.
5. EX-55 : For posterior tongue surface transverse / vertical position deep cut, it is easy to coronal the depth.
6. TR-33 : Main diamond burs, reduce the remaining tooth structure, prepare forming axial plane, shoulder, coronal edge and scotia side.
7. TR-D10 : Remove maxillofacial position depth, prepare the original maxillofacial form.
8. FO-23 : Trim lingual fossa, prepare origin lingual fossa shape.
9. SF-D41F : Trim shoulder and axial plane.
10. TR-08F : Trim shoulder, remove the edge.
11. TR-D36F : Trim axial plane, Maxillofacial, coronal edge and shoulder. This diamond bur's angle is large.
12. SH 3414 : Use slow bent nose to polish axial plane, Maxillofacial. (Use water to polish under the under vital pulp, to prevent the heat stimulate Tooth nerve.)


Kit of Preparation

Kit Of Preparation Teeth for Porcelain veneers

Order No: FG 0807D

Packing : 8 pcs/kit

Item No.: BR-12,EX-56,EX-57,EX-58,TR-32,SO-20,TR-D36F,SH 3414



Introduction


1. BR-12 : For the edge of the gingival to removing the depth of 0.3-0.5MM and extend to the collar near the distal and cutting edge.(Window)
2. EX-56 : Removing labial surface depth 0.3MM instructions ditch.
3. EX-57 : Removing labial surface depth 0.5MM instructions ditch.
4. EX-58 : The longest diamond bur is 1.5MM, for cut end located,but tand encasement mode. It is easy to control.
5. TR-32 : Main diamond burs, reduce the remaining tooth structure, prepare cut end removing 0.5-1.0MM (butt mode)removing 0.5-0.8MM shoulder for surrounding the lingual portion of the tooth.
6. SO-20 : The diamond burs extend to the collar near the distal, removing labial surface 0.5-0.7MM, removal undercuts and trimming interproximal. Please do not destruct interproximal contact point.
7. TR-D36F : Finishing every ends, removing sharp line angle for preparation and no angle shoulder.
8. SH 3414 : Use slow bent nose to polish axial plane, Maxillofacial. (Use water to polish under the under vital pulp, to prevent the heat stimulate tooth nerve.)

Kit Of Preparation Teeth for Porcelain veneers

Order No:FG 0811D

Packing:11 pcs/kit

Item No.:TR-33、WR-25、EX-60、EX-59、BR-12、TR-D33F、TR-D36F、TB1570、TR-31、SH 3414、TR-86F




Introduction

- 1.TR-33: For depth of cut end ; Removing labial surface.
- 2.WR-25: For removing cut end.
- 3.EX-60: For veneer depth 0.3mm.
- 4.EX-59: For veneer depth 0.5mm.
- 5.BR-12: For the Shoulder depth.
- 6.TR-D33F: For grinding the labial surface of the teeth.
- 7.TR-D36F: For Polishing the labial surface of the teeth.
- 8.TB1570: For forming the shoulder and finishing the labial surface.
- 9.TR-31: For the preparation of the shape like '<' between the adjacent faces .
- 10.SH3414: For Fine polishing of cut ends and labial surface.
- 11.TR-86F: Used for finishing and polishing the lingual edge of the teeth after bonding.

Kit of Preparation

Kit of Inlay/onlay preparation

Order No： FG2010D
Packing： 10pcs/kit
Item No: TF-08, TF-D15, TF-D15F, TF-06, TF-07, TF-07F, TR-21, SF-31, SF-31F, SH 3114m



Introduction

- 1.TF-08: It is used to develop the inlay preparation of the maxillofacial fossa and groove of the posterior teeth with a width of 2.0-2.2mm and a depth of 2.0-2.2mm (mainly premolars).
2. TF-D15： It is used to prepare the maxillofacial fossa and groove of the posterior teeth with a width of 2.2-5mm and depth and the inlay of the functional tip (mainly premolars).
3. TF-D15F： It is used for fine repair of 2.0-5mm wide and deep posterior maxillofacial fossa and groove and functional tip inlay preparation (mainly premolars).
- 4.TF-06: It is used to develop the inlay preparation of 2.0-2.4mm width and depth of the maxillofacial fossa and groove of the posterior teeth and the functional tip (mainly molars).
5. TF-07:It is used to develop the inlay preparation of 2.4-4mm width and depth of the posterior maxillofacial fossa and groove and functional tip (mainly molars).
6. TF-07F:It is used for fine repair of 2.4-4mm wide and deep maxillofacial fossa and groove of posterior teeth and inlay preparation of functional tip (mainly molars).
7. TR-21:It is used to open the cheek of the inlay, the axial wall of the lingual adjacent surface to the self-cleaning area.
8. SF-31:or preparing and shaping 1.0mm shoulder stand of the inlay gingival step area.
9. SF-31F:or fine trimming 1.0mm shoulder stand of the inlay gingival step area.
- 10.SH 3114m:It is used to polish the enamel edge of the prepared maxillofacial teeth.

Pulp Removal Burs Kit

Order No： FG5008D
Packing： 8pcs/kit
Item No: BR-13, BR-L26, BR-L21, BR-L25, EX-10, EX-30, SR-L01, DD3873



Introduction

1. BR-13:It is used for anterior teeth, removal of old and decayed teeth and pulp opening
2. BR-L26:It is used for molars, removal of old and decayed teeth and pulp opening
- 3.BR-L21: It is used for minimally invasive, removing old and decaying, and penetrating dental pulp cavity
4. BR-L25:Used to remove the crown and pulp
5. EX-10:The safety head is designed to be able to contact the bottom without damage, open 4 degrees outside, and trim the side wall
6. EX-30:The safety head is designed to be able to contact the bottom without damage, open 4 degrees outside, and trim the side wall
7. SR-L01:The safety head design can be placed at the root canal mouth, open to eliminate resistance, and smooth the side wall
8. DD3873:Crown pulp removal

Kit of Preparation

Grit : Mid-Coarse

Packing(pcs): 5

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 806 314...

G/BR-12	G/BR-15	G/WR-25	G/EX-49	G/EX-50	G/EX-52	G/EX-54	G/EX-55	G/EX-56	G/EX-57	G/EX-58	G/EX-59
$\emptyset 1.1$	$\emptyset 1.7$	3.1x0.9	1.7x8	2.2x8	1.3x2.2	3.1x4	3.1x3	1.7x6	2.2x6	1.1x1.6	2.1x7
001 524 011	001 524 017	068 524 031	522 524 017	522 524 022	500 524 013	552 524 031	552 524 031	552 524 017	552 524 022	500 524 011	522 524 021

Grit : Mid-Coarse

Packing(pcs): 5

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 806 314...

G/EX-60	G/FO-23	G/FO-41	G/SO-20	G/TC-34	G/TF-06	G/TF-07	G/TF-08	G/TF-D15	G/TR-03	G/TR-05	G/TR-09	G/TR-21
1.9x7	1.9x3.8	1.7x4	1.1x8	1.6x10	2.3x2.2	2.4x4	2x2.2	2x5	2.1x4	2.4x4	2.5x4	1.5x6
522 524 019	257 524 019	257 524 017	289 524 011	166 524 016	205 524 023	205 524 024	205 524 020	170 524 020	196 524 021	196 524 024	196 524 025	197 524 015

Grit : Mid-Coarse

Packing(pcs): 5

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 806 314...

G/TR-22	G/TR-31	G/TR-32	G/TR-33	G/TR-34	G/TR-86	G/TR-D10	G/TR-D36	G/SF-31	G/SF-D31	G/SF-45	G/SR-26	G/SR-41
1.6x6	1.2x8	1.5x8	1.6x8	1.8x8	2.5x5	2x5	2.5x7	1.2x6	1.3x5	1.4x8	1.3x8	1.3x6
197 524 016	198 524 012	198 524 015	198 524 016	198 524 018	277 524 025	196 524 020	197 524 025	110 524 012	110 524 013	108 514 014	141 524 013	139 524 013

Grit : Coarse

Packing(pcs): 5

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 806 314...

						Fine					
G/BR-31C	G/BR-32C	G/EX-84C	G/EX-85C	G/TR-39C	G/TR-D23C	G/BR-S11F	G/FO-S12F	G/FO-41F	G/SF-31F	G/SF-45F	G/SF-D41F
$\emptyset 1.6$	$\emptyset 1.3$	1.4x1.6	1.4x1.2	1.9x8	1.7x6	$\emptyset 0.7$	0.8x4	1.5x4	1.1x6	1.2x8	1.1x8
001 534 016	001 534 013	108 534 014	108 534 014	197 534 019	197 534 017	001 514 007	247 514 008	257 514 015	110 514 011	146 514 012	111 514 011

Grit : Fine

Packing(pcs): 5

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 806 314...

G/SR-S11F	G/TC-S11F	G/TC-09F	G/TF-07F	G/TF-D15F	G/TR-03F	G/TR-05F	G/TR-08F	G/TR-17F	G/TR-33F	G/TR-34F	G/TR-86F
0.8x2.2	0.7x2.2	0.6x5	2.2x4	1.9x5	2x4	2.3x4	1.4x1	2.1x5	1.6x8	1.7x8	2.3x5
138 514 008	161 514 007	164 514 006	205 514 022	170 514 019	196 514 020	196 514 023	537 514 014	196 514 021	198 514 016	198 514 017	277 514 023

Grit : Fine

Packing(pcs): 5

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 806 314...

		Extra Fine										
G/TR-D33F	G/TR-D36F			G/TC-21EF	G/TR-03EF	G/TR-05EF	G/TR-34EF	G/TR-D36EF	G/TR-86EF	G/EX-41EF	G/SF-45EF	G/FO-18EF
1.5x8	2.3x7			0.8x6	1.9x4	2.2x4	1.6x8	2.3x7	2.3x5	1x0.5	1.2x8	1.2x10
198 514 015	198 514 023			164 504 008	196 504 019	196 504 022	198 504 016	198 504 023	277 504 023	150 504 010	146 504 012	250 504 012

CD-3 Ceramic Diamond Grinders



Ceramic Diamond Grinders
Special for zirconia

Material: diamond grits mixing with ceramic powder.

Usage: For Zirconia inner crow fast finishing dry grinding without water cooling. (It is a effective fast way to instead of the diamond burs with water cooling way)

Caution: Please use the products in allowed RPM. OR it may Take the products damaged!

																		
Grit :	Mid-Coarse																	Fine
Packing(pcs):	6/12																	
Order No. :		CD 351 3	CD 324 3	CD 314 3	CD 312 3	CD 306 3	CD 301 3	CD 314 4	CD 306 4									
Recommended R.P.M.:		10,000~15,000 / min																
Size $\Phi \times L$ (mm):		$\Phi 3$	3.6 x 7	3 x 7	4 x 10.5	3 x 7	3 x 7	3 x 7	3 x 7									
ISO No. : 805...		204 001 524 030	204 237 524 036	204 197 524 030	204 174 524 040	204 273 524 030	204 110 524 030	204 197 514 030	204 273 514 030									

RD-1/3 Rubber Diamond Polishers



Rubber Diamond Polishers
Special for dental clinic porcelain and zirconia.

Material: diamond grits mixing with rubber.

Usage: For supper bright polishing for porcelain, without glaze.

Caution: Please use the products in allowed RPM. OR it may Take the products damaged!

									
Grit :	Mid-Coarse			Fine			Extra Fine		
Packing(pcs): 6/12/100									
Order No. :	RD 101 3	RD 111 3	RD 113 3	RD 101 4	RD 111 4	RD 113 4	RD 101 5	RD 111 5	RD 113 5
Recommended R.P.M.:	15,000~25,000 / min			15,000~20,000 / min			10,000~15,000 / min		
Size $\Phi \times L$ (mm):	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
ISO No. : 802...	314 110 524 030	314 243 524 030	314 243 524 050	314 110 514 030	314 243 514 030	314 243 514 050	314 110 504 030	314 243 504 030	314 243 504 050

Grit : Mid-Coarse							
Packing(pcs): 6/12/100							
Order No. :	RD 301 3	RD 311 3	RD 312 3	RD 313 3	RD 316 3	RD 320 3	RD 321 3
Recommended R.P.M.:	15,000~25,000 / min						
Size $\Phi \times L$ (mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5
ISO No. : 802...	204 110 524 030	204 243 524 030	204 243 524 040	204 243 524 050	204 174 524 035	204 030 524 060	204 030 524 090

RD-1/3 Rubber Diamond Polishers



Rubber Diamond Polishers
Special for dental clinic porcelain and zirconia.

Material: diamond grits mixing with rubber.

Usage: For supper bright polishing for porcelain, without glaze.

Caution: Please use the products in allowed RPM. OR it may Take the products damaged!

Grit : Mid-Coarse			
Packing(pcs): 6/12/100			
Order No. :	RD 331 3	RD 332 3	RD 341 3
Recommended R.P.M.:	15,000~25,000 / min		
Size $\Phi \times L$ (mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 802...	204 371 524 100	204 041 524 100	204 304 524 100

Grit : Fine							
Packing(pcs): 6/12/100							
Order No. :	RD 301 4	RD 311 4	RD 312 4	RD 313 4	RD 316 4	RD 320 4	RD 321 4
Recommended R.P.M.:	15,000~20,000 / min						
Size $\Phi \times L$ (mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5
ISO No. : 802...	204 110 514 030	204 243 514 030	204 243 514 040	204 243 514 050	204 174 514 035	204 030 514 060	204 030 514 090

Grit : Fine			
Packing(pcs): 6/12/100			
Order No. :	RD 331 4	RD 332 4	RD 341 4
Recommended R.P.M.:	15,000~20,000 / min		
Size $\Phi \times L$ (mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 802...	204 371 514 100	204 041 514 100	204 304 514 100

Grit : Extra Fine							
Packing(pcs): 6/12/100							
Order No. :	RD 301 5	RD 311 5	RD 312 5	RD 313 5	RD 316 5	RD 320 5	RD 321 5
Recommended R.P.M.:	10,000~15,000 / min						
Size $\Phi \times L$ (mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9 x 8.5
ISO No. : 802...	204 110 504 030	204 243 504 030	204 243 504 040	204 243 504 050	204 174 504 035	204 030 504 060	204 030 504 090

Grit : Extra Fine

Packing(pcs): 6/12

Order No. :

Recommended R.P.M.:

Size $\Phi \times L$ (mm):

ISO No. : 802...



RD 331 5



RD 332 5



RD 341 5

10,000~15,000 / min

10 x 1

10 x 2

10 x 1.8

204 371 504 100

204 041 504 100

204 304 504 100



SD-1/3

Silicone
Diamond Polishers

Silicone Diamond Polishers
Special for dental clinic resin.

Material:diamond grits mixing with silicone rubber.

Usage: For repairing and polishing resin teeth/hidden
artificial teeth.

Caution: Please use the products in allowed RPM.OR it
mayTake the products damaged!

Grit :

Extra Fine

Packing(pcs): 6/12/100

Order No. :

Recommended R.P.M.:

Size $\Phi \times L$ (mm):

ISO No. : 802...



SD 101 5



SD 111 5



SD 113 5



SD 101 6



SD 111 6



SD 113 6

10,000~12,000 / min

3 x 7

3 x 7

5 x 10

3 x 7

3 x 7

5 x 10

314 110 504 030

314 243 504 030

314 243 504 050

314 110 494 030

314 243 494 030

314 243 494 050

Ultra Fine

10,000~12,000 / min

Silicone Diamond Polishers
Special for dental clinic resin.

Material:diamond grits mixing with silicone rubber.

Usage: For repairing and polishing resin teeth/hidden
artificial teeth.

Caution: Please use the products in allowed RPM.OR it may
Take the products damaged!



SD-1/3

Silicone
Diamond Polishers



Grit : Extra Fine

Packing(pcs): 6/12

Order No. :

Recommended R.P.M.:

Size $\Phi \times L$ (mm):

ISO No. : 808...



SD 301 5



SD 311 5



SD 312 5



SD 313 5



SD 316 5



SD 320 5



SD 321 5

10,000~12,000 / min

3 x 7

3 x 7

4 x 12

5 x 10

3.5 x 10.5

6 x 9.5

9 x 8.5

204 110 504 030

204 243 504 030

204 243 504 040

204 243 504 050

204 174 504 035

204 030 504 060

204 030 504 090



SD 331 5



SD 332 5



SD 341 5

10,000~12,000 / min

10 x 1

10 x 2

10 x 1.8

204 371 504 100

204 041 504 100

204 304 504 100

Grit : Extra Fine

Packing(pcs): 6/12

Order No. :

Recommended R.P.M.:

Size $\Phi \times L$ (mm):

ISO No. : 808...



SD 301 6



SD 311 6



SD 312 6



SD 313 6



SD 316 6



SD 320 6



SD 321 6

10,000~12,000 / min

3 x 7

3 x 7

4 x 12

5 x 10

3.5 x 10.5

6 x 9.5

9 x 8.5

204 110 494 030

204 243 494 030

204 243 494 040

204 243 494 050

204 174 494 035

204 030 494 060

204 030 494 090



SD 331 6



SD 332 6



SD 341 6

10,000~12,000 / min

10 x 1

10 x 2

10 x 1.8

204 371 494 100

204 041 494 100

204 304 494 100

Grit : Ultra Fine

Packing(pcs): 6/12

Order No. :

Recommended R.P.M.:

Size $\Phi \times L$ (mm):

ISO No. : 808...

RF-1/3
Rubber Polishers

Rubber Polishers
Special for dental clinic Non-precious and semi-Precious metal
Material:rubber mixing with normal abrasive.
Usage:For repairing dental clinic silver-mercury alloy/titanium alloy/nickel-chrome/cobalt-chrome alloy.
Caution: Please use the products in allowed RPM.OR it may take the products damaged!



Grit :	Coarse			Fine		
Packing(pcs): 6/12/100						
Order No. :	RF 101 2	RF 111 2	RF 113 2	RF 101 4	RF 111 4	RF 113 4
Recommended R.P.M:	10,000~15,000 / min			10,000~15,000 / min		
Size Φ × L(mm):	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
ISO No. : 652...	314 110 534 030	314 243 534 030	314 243 534 050	314 110 514 030	314 243 514 030	314 243 514 050

Grit :	Extra Fine			Ultra Fine		
Packing(pcs): 6/12/100						
Order No. :	RF 101 5	RF 111 5	RF 113 5	RF 101 6	RF 111 6	RF 113 6
Recommended R.P.M:	10,000~15,000 / min			10,000~15,000 / min		
Size Φ × L(mm):	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
ISO No. : 652...	314 110 504 030	314 243 504 030	314 243 504 050	314 110 494 030	314 243 494 030	314 243 494 050

Grit : Coarse							
Packing(pcs): 6/12/100							
Order No. :	RF 301 2	RF 311 2	RF 312 2	RF 313 2	RF 316 2	RF 320 2	RF 321 2
Recommended R.P.M:	10,000~15,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 652...	204 110 534 030	204 243 534 030	204 243 534 040	204 243 534 050	204 174 534 035	204 030 534 060	204 030 534 090

Grit : Coarse			
Packing(pcs): 6/12/100			
Order No. :	RF 331 2	RF 332 2	RF 341 2
Recommended R.P.M:	10,000~15,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 652...	204 371 534 100	204 041 534 100	204 304 534 100

Grit : Fine							
Packing(pcs): 6/12/100							
Order No. :	RF 301 4	RF 311 4	RF 312 4	RF 313 4	RF 316 4	RF 320 4	RF 321 4
Recommended R.P.M:	10,000~15,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 652...	204 110 514 030	204 243 514 030	204 243 534 040	204 243 514 050	204 174 514 035	204 030 514 060	204 030 514 090

Rubber Polishers
Special for dental clinic Non-precious and semi-Precious metal
Material:rubber mixing with normal abrasive.
Usage:For repairing dental clinic silver-mercury alloy/titanium alloy/nickel-chrome/cobalt-chrome alloy.
Caution: Please use the products in allowed RPM.OR it may Take the products damaged!



RF-3
Rubber Polishers

Grit : Fine			
Packing(pcs): 6/12/100			
Order No. :	RF 331 4	RF 332 4	RF 341 4
Recommended R.P.M:	10,000~15,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 652...	204 371 514 100	204 041 514 100	204 304 514 100

Grit : Extra Fine							
Packing(pcs): 6/12/100							
Order No. :	RF 301 5	RF 311 5	RF 312 5	RF 313 5	RF 316 5	RF 320 5	RF 321 5
Recommended R.P.M:	10,000~15,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 652...	204 110 504 030	204 243 504 030	204 243 504 040	204 243 504 050	204 174 504 035	204 030 504 060	204 030 504 090

Grit : Extra Fine			
Packing(pcs): 6/12/100			
Order No. :	RF 331 5	RF 332 5	RF 341 5
Recommended R.P.M:	10,000~15,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 652...	204 371 504 100	204 041 504 100	204 304 504 100

Grit : Ultra Fine							
Packing(pcs): 6/12/100							
Order No. :	RF 301 6	RF 311 6	RF 312 6	RF 313 6	RF 316 6	RF 320 6	RF 321 6
Recommended R.P.M:	10,000~12,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 652...	204 110 494 030	204 243 494 030	204 243 494 040	204 243 494 050	204 174 494 035	204 030 494 060	204 030 494 090

Grit : Ultra Fine			
Packing(pcs): 6/12/100			
Order No. :	RF 331 6	RF 332 6	RF 341 6
Recommended R.P.M:	10,000~12,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 652...	204 371 494 100	204 041 494 100	204 304 494 100

SH-1

Silicone Polishers



Silicone Polishers
Special for porcelain/resin teeth.

Material:silicon rubber mixing with normal abrasive.

Usage:For smoothing and initial polishing porcelain/resin teeth.
(super bright or without glaze polishing suggest use RD2/0 and SD2/0 series product)

Caution: Please use the products in allowed RPM.OR it may
Take the products damaged!

Grit :	Coarse			Mid-Coarse			Fine		
Packing(pcs): 6/12/100									
Order No. :	SH 101 2	SH 111 2	SH 113 2	SH 101 3	SH 111 3	SH 113 3	SH 101 4	SH 111 4	SH 113 4
Recommended R.P.M:	10,000~15,000 / min			10,000~15,000 / min			10,000~15,000 / min		
Size Φ × L(mm):	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
ISO No. : 658...	314 110 534 030	314 243 534 030	314 243 534 050	314 110 524 030	314 243 524 030	314 243 524 050	314 110 514 030	314 243 514 030	314 243 514 050

SH-3

Silicone Polishers

							
Grit : Mid-Coarse	SH 301 6B	SH 311 6B	SH 312 6B	SH 313 6B	SH 316 6B	SH 320 6B	SH 321 6B
Packing(pcs): 6/12/100							
Order No. :							
Recommended R.P.M:				10,000~12,000 / min			
Size Φ x L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 524 030	204 243 524 030	204 243 524 040	204 243 524 050	204 174 524 035	204 030 524 060	204 030 524 090

			
Grit : Mid-Coarse			
Packing(pcs): 6/12/100			
Order No. :	SH 331 6B	SH 332 6B	SH 341 6B
Recommended R.P.M:	10,000~12,000 / min		
Size Φ x L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 658...	204 371 524 100	204 041 524 100	204 304 524 100

SH-3

Silicone Polishers



Silicone Polishers
Special for porcelain/resin teeth.

Material:silicon rubber mixing with normal abrasive.

Usage:For smoothing and initial polishing porcelain/resin teeth.
(super bright or without glaze polishing suggest use RD2/0 and SD2/0 series product)

Caution: Please use the products in allowed RPM.OR it may
Take the products damaged!

							
Grit : Coarse	SH 301 2	SH 311 2	SH 312 2	SH 313 2	SH 316 2	SH 320 2	SH 321 2
Packing(pcs): 6/12/100							
Order No. :							
Recommended R.P.M:				10,000~15,000 / min			
Size Φ x L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 534 030	204 243 534 030	204 243 534 040	204 243 534 050	204 174 534 035	204 030 534 060	204 030 534 090

			
Grit : Coarse			
Packing(pcs): 6/12/100			
Order No. :	SH 331 2	SH 332 2	SH 341 2
Recommended R.P.M:	10,000~15,000 / min		
Size Φ x L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 658...	204 371 534 100	204 041 534 100	204 304 534 100

							
Grit : Mid-Coarse							
Packing(pcs): 6/12/100							
Order No. :	SH 301 3	SH 311 3	SH 312 3	SH 313 3	SH 316 3	SH 320 3	SH 321 3
Recommended R.P.M:				10,000~12,000 / min			
Size Φ x L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 524 030	204 243 524 030	204 243 524 040	204 243 524 050	204 174 524 035	204 030 524 060	204 030 524 090

			
Grit : Mid-Coarse	SH 331 3	SH 332 3	SH 341 3
Packing(pcs): 6/12/100			
Order No. :			
Recommended R.P.M:		10,000~12,000 / min	
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 658...	204 371 524 100	204 041 524 100	204 304 524 100

SH-3

Silicone Polishers

Silicone Polishers
Special for porcelain/resin teeth.

Material:silicon rubber mixing with normal abrasive.

Usage:For smoothing and initial polishing porcelain/resin teeth.
(super bright or without glaze polishing suggest use RD2/0 and SD2/0 series product)

Caution: Please use the products in allowed RPM.OR it may
Take the products damaged!

Grit : Fine							
Packing(pcs): 6/12/100							
Order No. :	SH 301 4D	SH 311 4D	SH 312 4D	SH 313 4D	SH 316 4D	SH 320 4D	SH 321 4D
Recommended R.P.M:	10,000~12,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 514 030	204 243 514 030	204 243 514 040	204 243 514 050	204 174 514 035	204 030 514 060	204 030 514 090

Grit : Fine			
Packing(pcs): 6/12/100			
Order No. :	SH 331 4D	SH 332 4D	SH 341 4D
Recommended R.P.M:	10,000~12,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 658...	204 371 514 100	204 041 514 100	204 304 514 100

Grit : Fine							
Packing(pcs): 6/12/100							
Order No. :	SH 301 4	SH 311 4	SH 312 4	SH 313 4	SH 316 4	SH 320 4	SH 321 4
Recommended R.P.M:	10,000~12,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 514 030	204 243 514 030	204 243 514 040	204 243 514 050	204 174 514 035	204 030 514 060	204 030 514 090

Grit : Fine			
Packing(pcs): 6/12/100			
Order No. :	SH 331 4	SH 332 4	SH 341 4
Recommended R.P.M:	10,000~12,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 658...	204 371 514 100	204 041 514 100	204 304 514 100

Silicone Polishers
Special for precious metal workpiece.

Material:silicon rubber mixing with normal abrasive.

Usage:For finishing and polishing precious/semi-precious dental workpiece,such as finishing,smoothing and polishing golden teeth and golden inner crown.

Caution: Please use the products in allowed RPM.OR it may
Take the products damaged!

Grit :									
Packing(pcs): 6/12/100									
Order No. :	SG 101 2	SG 111 2	SG 113 2	SG 101 3	SG 111 3	SG 113 3	SG 101 4	SG 111 4	SG 113 4
Recommended R.P.M:	7,000~12,000 / min			7,000~12,000 / min			7,000~12,000 / min		
Size Φ × L(mm):	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10	3 x 7	3 x 7	5 x 10
ISO No. : 658...	314 110 534 030	314 243 534 030	314 243 534 050	314 110 524 030	314 243 524 030	314 243 524 050	314 110 514 030	314 243 514 030	314 243 514 050

Grit : Coarse							
Packing(pcs): 6/12/100							
Order No. :	SG 301 2	SG 311 2	SG 312 2	SG 313 2	SG 316 2	SG 320 2	SG 321 2
Recommended R.P.M:	10,000~15,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 534 030	204 243 534 030	204 243 534 040	204 243 534 050	204 174 534 035	204 030 534 060	204 030 534 090

Grit : Coarse							
Packing(pcs): 6/12/100							
Order No. :	SG 301 2	SG 311 2	SG 312 2	SG 313 2	SG 316 2	SG 320 2	SG 321 2
Recommended R.P.M:	10,000~15,000 / min						
Size Φ × L(mm):							
ISO No. : 658...							

SG-1

Silicone Polishers



SG-3

Silicone Polishers

SG-3
Silicone Polishers

Silicone Polishers
Special for precious metal workpiece.
Material:silicon rubber mixing with normal abrasive.
Usage:For finishing and polishing precious/semi-precious dental workpiece,such as finishing,smoothing and polishing golden teeth and golden inner crown.
Caution: Please use the products in allowed RPM.OR it may Take the products damaged!

Grit : Mid-Coarse							
Packing(pcs): 6/12/100							
Order No. :							
Recommended R.P.M:	10,000~15,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 524 030	204 243 524 030	204 243 524 040	204 243 524 050	204 174 524 035	204 030 524 060	204 030 524 090

Grit : Mid-Coarse			
Packing(pcs): 6/12/100			
Order No. :			
Recommended R.P.M:	10,000~15,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 658...	204 371 524 100	204 041 524 100	204 304 524 100

Grit : Fine							
Packing(pcs): 6/12/100							
Order No. :							
Recommended R.P.M:	7,000~12,000 / min						
Size Φ × L(mm):	3 x 7	3 x 7	4 x 12	5 x 10	3.5 x 10.5	6 x 9.5	9.5 x 8
ISO No. : 658...	204 110 514 030	204 243 514 030	204 243 514 040	204 243 514 050	204 174 514 035	204 030 514 060	204 030 514 090

Grit : Fine			
Packing(pcs): 6/12/100			
Order No. :			
Recommended R.P.M:	7,000~12,000 / min		
Size Φ × L(mm):	10 x 1	10 x 2	10 x 1.8
ISO No. : 658...	204 371 514 100	204 041 514 100	204 304 514 100

CG-1/3 CW-1/3
Ceramic Grinders

Ceramic Grinder
Special for dental clinic metal,resin and porcelain.
Material:silicon carbide mixing with cermic powder.
Usage:For grinding dental clinic used titanium alloy/pure titanium/nickel-cobalt alloy/resin.
Caution: Please use the products in allowed RPM.OR it may Take the products damaged!

Grit : Coarse				
Packing(pcs): 6/12/100				
Order No. :				
Recommended R.P.M:	7,000~10,000 / min			
Size Φ × L(mm):	3 x 3	2.5 x 6	2.8 x 7	3 x 7
ISO No. : 655...	314 001 533 030	314 171 533 025	314 161 533 028	314 273 533 030

Grit : Coarse								
Packing(pcs): 6/12/100								
Order No. :								
Recommended R.P.M:	7,000~10,000 / min							
Size Φ × L(mm):	3 x 3	8 x 2	6 x 1.5	3.6 x 7	2.5 x 6	2.8 x 7	3 x 7	2.5 x 6
ISO No. : 655...	204 001 533 030	204 304 533 080	204 041 533 060	204 237 533 036	204 171 533 025	204 161 533 028	204 273 533 030	204 110 533 025

Grit : Fine				
Packing(pcs): 6/12/100				
Order No. :				
Recommended R.P.M:	7,000~10,000 / min			
Size Φ × L(mm):	3 x 3	2.5 x 6	2.8 x 7	3 x 7
ISO No. : 635...	314 001 513 030	314 171 513 025	314 161 513 028	314 273 513 030

Grit : Fine								
Packing(pcs): 6/12/100								
Order No. :								
Recommended R.P.M:	7,000~10,000 / min							
Size Φ × L(mm):	3 x 3	8 x 2	6 x 1.5	3.6 x 7	2.5 x 6	2.8 x 7	3 x 7	2.5 x 6
ISO No. : 635...	204 001 513 030	204 304 513 080	204 041 513 060	204 237 513 036	204 171 513 025	204 161 513 028	204 273 513 030	204 110 513 025

CD-0/2
Ceramic
Diamond Grinders

Ceramic Diamond Grinders
Special for zirconia

Material:diamond grits mixing with ceramic powder.

Usage:For zirconia inner crown fast finishing, dry grinding,
without water cooling (it is a effective fast way to instead
Of the diamond burs with water cooling way)

Caution: Please use the products in allowed RPM. OR it
mayTake the products damaged!

Rubber Diamond Polishers
Special for Zirconia/porcelain workpiece

Material:Diamond grits mixing with synthetic rubber

Usage:For smoothing and polishing zirconia inner crown,
zirconia veneer/porcelain teeth fast super bright .(Avoid the
color differents and porcelain fragile caused by second glaze
and baking.And it also can save time)

Caution: Please use the products in allowed RPM. OR it
mayTake the products damaged!

RD-0
Rubber
Diamond Polishers

Grit : Mid-Coarse

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 805...

						
CD 001 3	CD 006 3	CD 008 3	CD 007 3	CD 009 3	CD 011 3	CD 012 3
15,000-30,000 / min						
13 x 7	16 x 2	22 x 2	16 x 3.5	22 x 4.5	16 x 1.5	22 x 2
900 030 524 130	900 372 524 160	900 372 524 220	900 041 524 160	900 041 524 220	900 307 524 130	900 307 524 220

Grit : Mid-Coarse

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 805...

						
CD 201 3	CD 202 3	CD 203 3	CD 212 3	CD 214 3	CD 221 3	CD 222 3
20,000-35,000 / min						
3 x 7	5 x 12	7 x 12	4 x 10.5	3 x 7	5 x 2.5	13 x 3.5
104 110 524 030	104 113 524 050	104 113 524 070	104 174 524 040	104 197 524 030	104 014 524 050	104 024 524 130

Mid-Coarse

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 805...

				
CD 224 3	CD 225 3	CD 231 3	CD 232 3	CD 251 3
20,000-35,000 / min				
3.6 x 7	6 x 8	12 x 2	12 x 3	3 x 3
104 237 524 036	104 030 524 060	104 041 524 120	104 042 524 120	104 001 524 030

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 805...

			
CD 232 4	CD 202 4	CD 212 4	CD 214 4
20,000-35,000 / min			
12 x 3	5 x 12	4 x 10.5	3 x 7
104 042 514 120	104 113 514 050	104 174 514 040	104 174 514 030

Grit :

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 802...

				
RD 011 1	RD 011 3	RD 011 4	RD 011 5	RD 011 6
15,000-25,000 / min				
19 x 2	19 x 2	19 x 2	19 x 2	19 x 2
900 307 544 190	900 307 524 190	900 307 514 190	900 307 504 190	900 307 494 190

Grit :

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 802...

				
RD 012 1	RD 012 3	RD 012 4	RD 012 5	RD 012 6
15,000-25,000 / min				
26 x 3	26 x 3	26 x 3	26 x 3	26 x 3
900 307 544 260	900 307 524 260	900 307 514 260	900 307 504 260	900 307 494 260

Grit : Mid-Coarse

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 802...

						
RD 201 3	RD 202 3	RD 211 3	RD 212 3	RD 216 3	RD 221 3	RD 222 3
15,000-25,000 / min						
3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	18 x 3
104 110 524 030	104 113 524 050	104 243 524 030	104 243 524 040	104 174 524 040	104 030 524 090	104 026 524 180

Grit : Mid-Coarse

Packing(pcs): 1

Order No. :

Recommended R.P.M:

Size $\Phi \times L$ (mm):

ISO No. : 802...

				
RD 231 3	RD 232 3	RD 233 3	RD 241 3	RD 237 3A
15,000-25,000 / min				
12 x 1	12 x 2	12 x 4	12 x 2	15 x 2.5
104 371 524 120	104 041 524 120	104 042 524 120	104 304 524 120	104 041 524 150

RD-2
Rubber
Diamond Polishers

Rubber Diamond Polishers
Special for Zirconia/porcelain workpiece
Material:Diamond grits mixing with synthetic rubber
Usage:For smoothing and polishing zirconia inner crown, zirconia veneer/porcelain teeth fast super bright .(Avoid the color differents and porcelain fragile caused by second glaze and baking.And it also can save time)
Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Rubber Diamond Polishers
Special for Zirconia/porcelain workpiece
Material:Diamond grits mixing with synthetic rubber
Usage:For smoothing and polishing zirconia inner crown, zirconia veneer/porcelain teeth fast super bright .(Avoid the color differents and porcelain fragile caused by second glaze and baking.And it also can save time)
Caution: Please use the products in allowed RPM.OR it may Take the products damaged!

RD-2
Rubber
Diamond Polishers

Grit :							
Packing(pcs): 1							
Order No. :							
Recommended R.P.M:	15,000~25,000 / min						
Size $\Phi \times L(mm)$:	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	18 x 3
ISO No. : 802...	104 110 514 030	104 113 514 050	104 243 514 030	104 243 514 040	104 174 514 040	104 030 514 090	104 026 514 180

Packing(pcs): 1					
Order No. :					
Recommended R.P.M:	15,000~25,000 / min				
Size $\Phi \times L(mm)$:	12 x 1	12 x 2	12 x 4	12 x 2	15x2.5
ISO No. : 802...	104 371 514 120	104 041 514 120	104 042 514 120	104 304 514 120	104 041 514 150

Extra Fine							
Packing(pcs): 1							
Order No. :							
Recommended R.P.M:	10,000~20,000 / min						
Size $\Phi \times L(mm)$:	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	18 x 3
ISO No. : 802...	104 110 504 030	104 113 504 050	104 243 504 030	104 243 504 040	104 174 504 040	104 030 504 090	104 026 504 180

Extra Fine				
Packing(pcs): 1				
Order No. :				
Recommended R.P.M:	10,000~20,000 / min			
Size $\Phi \times L(mm)$:	12 x 1	12 x 2	12 x 4	12 x 2
ISO No. : 802...	104 371 504 120	104 041 504 120	104 042 504 120	104 304 504 120

Grit : Ultra Fine							
Packing(pcs): 1							
Order No. :							
Recommended R.P.M:	10,000~20,000 / min						
Size $\Phi \times L(mm)$:	3 x 7	5 x 12	3 x 7	4 x 12	4 x 12	9 x 8.5	18 x 3
ISO No. : 802...	104 110 494 030	104 113 494 050	104 243 494 030	104 243 494 040	104 174 494 040	104 030 494 090	104 026 494 180

Grit : Ultra Fine				
Packing(pcs): 1				
Order No. :				
Recommended R.P.M:	10,000~20,000 / min			
Size $\Phi \times L(mm)$:	12 x 1	12 x 2	12 x 4	12 x 2
ISO No. : 802...	104 371 494 120	104 041 494 120	104 042 494 120	104 304 494 120

粒度 Grit : 超细 Ultra Fine				
包装 Packing(pcs): 1/12				
订货编号 Order No. :				
推荐转速 Recommended R.P.M:	15,000~25,000 / min			
规格 Size $\Phi \times L(mm)$:	4 x 12	6 x 16	15x2.5	15x2.5
ISO No. : 802...	104 243 494 040	104 243 494 060	104 041 494 150	104 041 494 150

SD-0/2
Silicone
Diamond Polishers

Extra Fine						
Packing(pcs): 1						
Order No. :						
Recommended R.P.M:	8,000~10,000 / min					
Size $\Phi \times L(mm)$:	22 x 3	12 x 1	4 x 12	22 x 3	12 x 1	4 x 12
ISO No. : 808...	900 372 504 220	104 371 504 120	104 243 504 040	900 372 494 220	104 371 494 120	104 243 494 040

CG-0/2 CW-0/2
Ceramic Grinders



Ceramic Grinder
Material:silicon carbide mixing with cermic powder.

Usage:For grinding titanium alloy/pure titanium/nickel-cobalt alloy .

Caution: Please use the products in allowed RPM.OR it may Take the products damaged!

Grit : Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M.:	15,000-30,000 / min						
Size $\Phi \times L$ (mm):	16 x 3	22 x 3	3 x 7	5 x 12	6.5 x 13	5 x 16	3 x 7
ISO No. : 655...	900 372 533 160	900 372 533 220	104 110 533 030	104 113 533 050	104 113 533 065	104 292 533 050	104 273 533 030

Grit : Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M.:	15,000-30,000 / min						
Size $\Phi \times L$ (mm):	4.5 x 7.5	2.8 x 7	2.5 x 6	3.5 x 10.5	5 x 10	5 x 2.5	7 x 3
ISO No. : 655...	104 257 533 045	104 161 533 028	104 171 533 025	104 173 533 035	104 266 533 050	104 014 533 050	104 014 533 070

Grit : Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M.:	15,000-25,000 / min						
Size $\Phi \times L$ (mm):	12.5 x 3.5	3.6 x 7	5 x 3	6 x 1.5	6.5 x 4	13 x 1.5	6.5 x 3
ISO No. : 655...	104 024 533 125	104 237 533 036	104 042 533 050	104 041 533 060	104 043 533 065	104 041 533 130	104 072 533 065

Grit : Coarse				
Packing(pcs): 12/100				
Order No. :				
Recommended R.P.M.:	15,000-25,000 / min			
Size $\Phi \times L$ (mm):	6.5 x 2	8 x 2	3 x 3	5 x 5
ISO No. : 655...	104 067 533 065	104 304 533 080	104 001 533 030	104 001 533 050

CG-0/2 CW-0/2
Ceramic Grinders



Ceramic Grinder
Material:silicon carbide mixing with cermic powder.

Usage:For grinding titanium alloy/pure titanium/nickel-cobalt alloy .

Caution: Please use the products in allowed RPM.OR it may Take the products damaged!

Grit : Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M.:	15,000-30,000 / min						
Size $\Phi \times L$ (mm):	16 x 3	22 x 3	3 x 7	5 x 12	6.5 x 13	5 x 16	3 x 7
ISO No. : 655...	900 372 533 160	900 372 533 220	104 110 533 030	104 113 533 050	104 113 533 065	104 292 533 050	104 273 533 030

Grit : Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M.:	15,000-30,000 / min						
Size $\Phi \times L$ (mm):	4.5 x 7.5	2.8 x 7	2.5 x 6	3.5 x 10.5	5 x 10	5 x 2.5	7 x 3
ISO No. : 655...	104 257 533 045	104 161 533 028	104 171 533 025	104 173 533 035	104 266 533 050	104 014 533 050	104 014 533 070

Grit : Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M.:	15,000-30,000 / min						
Size $\Phi \times L$ (mm):	12.5 x 3.5	3.6 x 7	5 x 3	6 x 1.5	6.5 x 4	13 x 1.5	6.5 x 3
ISO No. : 655...	104 024 533 125	104 237 533 036	104 042 533 050	104 041 533 060	104 043 533 065	104 041 533 130	104 072 533 065

Grit : Coarse				
Packing(pcs): 12/100				
Order No. :				
Recommended R.P.M.:	15,000-30,000 / min			
Size $\Phi \times L$ (mm):	6.5 x 2	8 x 2	3 x 3	5 x 5
ISO No. : 655...	104 067 533 065	104 304 533 080	104 001 533 030	104 001 533 050

Rubber
Polishers



Rubber Polishers
Special for Non-precious and semi-precious metals workpiece.
Material:synthetic rubber mixing with normalabrasive.
Usage:For finishing and polishing titanium alloy/pure titanium/nickel-chrome alloy/cobalt-chrome alloy dental workpiece,such as forsteel teeth/inner crown/steel frame /steel base finishing,smoothing and polishing.

Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Grit : Extra Coarse								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	20,000~30,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 652...	900 371 544 160	900 372 544 160	900 303 544 160	900 371 544 220	900 372 544 220	900 303 544 220	900 114 544 060	900 114 544 070

Grit : Coarse								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	20,000~30,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 652...	900 371 534 160	900 372 534 160	900 303 534 160	900 371 534 220	900 372 534 220	900 303 534 220	900 114 534 060	900 114 534 070

Grit : Mid-Coarse								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	20,000~30,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 652...	900 371 524 160	900 372 524 160	900 303 524 160	900 371 524 220	900 372 524 220	900 303 524 220	900 114 524 060	900 114 524 070

Grit :								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	10,000~20,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 652...	900 371 514 160	900 372 514 160	900 303 514 160	900 371 514 220	900 372 514 220	900 303 514 220	900 114 514 060	900 114 514 070

RF-2
Rubber
Polishers



Rubber Polishers
Special for Non-precious and semi-precious metals workpiece.
Material:synthetic rubber mixing with normalabrasive.
Usage:For finishing and polishing titanium alloy/pure titanium/nickel-chrome alloy/cobalt-chrome alloy dental workpiece,such as forsteel teeth/inner crown/steel frame /steel base finishing,smoothing and polishing.

Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Grit : Coarse									
Packing(pcs): 12/100									
Order No. :									
Recommended R.P.M:	20,000~35,000 / min								
Size Φ × L(mm):	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
ISO No. : 652...	104 110 534 030	104 113 534 050	104 292 534 050	104 243 534 030	104 243 534 040	104 243 534 060	104 031 534 060	104 030 534 090	104 026 534 180

Grit : Coarse									
Packing(pcs): 12/100									
Order No. :									
Recommended R.P.M:	20,000~35,000 / min								
Size Φ × L(mm):	12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2	
ISO No. : 652...	104 371 534 120	104 371 534 180	104 371 534 230	104 041 534 120	104 041 534 150	104 042 534 120	104 304 534 120	104 303 534 150	

Grit : Fine									
Packing(pcs): 12/100									
Order No. :									
Recommended R.P.M:	10,000~20,000 / min								
Size Φ × L(mm):	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
ISO No. : 652...	104 110 514 030	104 113 514 050	104 292 514 050	104 243 514 030	104 243 514 040	104 243 514 060	104 031 514 060	104 030 514 090	104 026 514 180

Grit : Fine									
Packing(pcs): 12/100									
Order No. :									
Recommended R.P.M:	10,000~20,000 / min								
Size Φ × L(mm):	12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2	
ISO No. : 652...	104 371 514 120	104 371 514 180	104 371 514 230	104 041 514 120	104 041 514 150	104 042 514 120	104 304 514 120	104 303 514 150	

RF-2
Rubber
Polishers



Rubber Polishers
Special for Non-precious and semi-precious metals workpiece.
Material:synthetic rubber mixing with normalabrasive.
Usage:For finishing and polishing titanium alloy/pure titanium/nickel-chrome alloy/cobalt-chrome alloy dental workpiece,such as forsteel teeth/inner crown/steel frame /steel base finishing,smoothing and polishing.

Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Grit : Extra Fine
Packing(pcs): 12/100
Order No. :
Recommended R.P.M:
Size $\Phi \times L$ (mm):
ISO No. : 652...

RF 201 5	RF 202 5	RF 206 5	RF 211 5	RF 212 5	RF 213 5	RF 216 5	RF 221 5	RF 222 5
10,000~12,000 / min								
3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
104 110 504 030	104 113 504 050	104 292 504 050	104 243 504 030	104 243 504 040	104 243 504 060	104 031 504 060	104 030 504 090	104 026 504 180

Grit : Extra Fine
Packing(pcs): 12/100
Order No. :
Recommended R.P.M:
Size $\Phi \times L$ (mm):
ISO No. : 652...

RF 231 5	RF 235 5	RF 236 5	RF 232 5	RF 234 5	RF 233 5	RF 241 5	RF 242 5
10,000~12,000 / min							
12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
104 371 504 120	104 371 504 180	104 371 504 230	104 041 504 120	104 041 504 150	104 042 504 120	104 304 504 120	104 303 504 150

Grit : Ultra Fine
Packing(pcs): 12/100
Order No. :
Recommended R.P.M:
Size $\Phi \times L$ (mm):
ISO No. : 652...

RF 201 6	RF 202 6	RF 206 6	RF 211 6	RF 212 6	RF 213 6	RF 216 6	RF 221 6	RF 222 6
10,000~12,000 / min								
3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	6 x 16	9 x 8.5	18 x 3
104 110 494 030	104 113 494 050	104 292 494 050	104 243 494 030	104 243 494 040	104 243 494 060	104 031 494 060	104 030 494 090	104 026 494 180

Grit : Ultra Fine
Packing(pcs): 12/100
Order No. :
Recommended R.P.M:
Size $\Phi \times L$ (mm):
ISO No. : 652...

RF 231 6	RF 235 6	RF 236 6	RF 232 6	RF 234 6	RF 233 6	RF 241 6	RF 242 6
10,000~12,000 / min							
12 x 1	18 x 1	23 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
104 371 494 120	104 371 494 180	104 371 494 230	104 041 494 120	104 041 494 150	104 042 494 120	104 304 494 120	104 303 494 150

Silicone Polishers
Special for porcelain/resin teeth.
Material:silicon rubber mixing with normal abrasive.
Usage:For smoothing and inital polishing porcelain/resin teeth.(super bright or without glaze polishing suggest use RD2/0 and SD2/0 series product)

Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Grit : Coarse
Packing(pcs): 100
Order No. :
Recommended R.P.M:
Size $\Phi \times L$ (mm):
ISO No. : 658...

SH 006 2	SH 007 2	SH 011 2	SH 008 2	SH 009 2	SH 012 2	SH 001 2	SH 002 2
15,000-30,000 / min							
16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
900 371 534 160	900 372 534 160	900 303 534 160	900 371 534 220	900 372 54 220	900 303 534 220	900 114 534 060	900 114 534 070

Grit : Mid-Coarse
Packing(pcs): 100
Order No. :
Recommended R.P.M:
Size $\Phi \times L$ (mm):
ISO No. : 658...

SH 006 3	SH 007 3	SH 011 3	SH 008 3	SH 009 3	SH 012 3	SH 001 3	SH 002 3
10,000~12,000 / min							
16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
900 371 524 160	900 372 524 160	900 303 524 160	900 371 524 220	900 372 524 220	900 303 524 220	900 114 524 060	900 114 524 070

Grit : Fine
Packing(pcs): 100
Order No. :
Recommended R.P.M:
Size $\Phi \times L$ (mm):
ISO No. : 658...

SH 006 4	SH 007 4	SH 011 4	SH 008 4	SH 009 4	SH 012 4	SH 001 4	SH 002 4
10,000~12,000 / min							
16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
900 371 514 160	900 372 514 160	900 303 514 160	900 371 514 220	900 372 514 220	900 303 514 220	900 114 514 060	900 114 514 070



SH-0
Silicone
Polishers

SH-2
Silicone
Polishers



Silicone Polishers
Special for porcelain/resin teeth.
Material:silicon rubber mixing with normal abrasive.
Usage:For smoothing and inital polishing porcelain/resin teeth.(super bright or without glaze polishing suggest use RD2/0 and SD2/0 series product)
Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Silicone Polishers
Special for porcelain/resin teeth.
Material:silicon rubber mixing with normal abrasive.
Usage:For smoothing and inital polishing porcelain/resin teeth.(super bright or without glaze polishing suggest use RD2/0 and SD2/0 series product)
Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

SH-2
Silicone
Polishers

Grit : Coarse									
Packing(pcs): 12/100									
Order No. :									
Recommended R.P.M:	15,000-30,000 / min								
Size Φ × L(mm):	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
ISO No. : 658...	104 110 534 030	104 113 534 050	104 292 534 050	104 243 534 030	104 243 534 040	104 243 534 060	104 174 534 040	104 030 534 090	104 026 534 180


Grit : Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M:	15,000-30,000 / min						
Size Φ × L(mm):	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
ISO No. : 658...	104 371 534 120	104 371 534 180	104 041 534 120	104 041 534 150	104 042 534 120	104 304 534 120	104 303 534 150

Grit : Mid-Coarse									
Packing(pcs): 12/100									
Order No. :									
Recommended R.P.M:	10,000-20,000 / min								
Size Φ × L(mm):	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
ISO No. : 658...	104 110 524 030	104 113 524 050	104 292 524 050	104 243 524 030	104 243 524 040	104 243 524 060	104 174 524 040	104 030 524 090	104 026 524 180

Grit : Mid-Coarse							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M:	10,000-20,000 / min						
Size Φ × L(mm):	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
ISO No. : 658...	104 371 524 120	104 371 524 180	104 041 524 120	104 041 524 150	104 042 524 120	104 304 524 120	104 303 524 150

Grit : Fine									
Packing(pcs): 12/100									
Order No. :									
Recommended R.P.M:	10,000-12,000 / min								
Size Φ × L(mm):	3 x 7	5 x 12	5 x 16	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
ISO No. : 658...	104 110 514 030	104 113 514 050	104 292 514 050	104 243 514 030	104 243 514 040	104 243 514 060	104 174 514 040	104 030 514 090	104 026 514 180

Grit : Fine							
Packing(pcs): 12/100							
Order No. :							
Recommended R.P.M:	10,000-12,000 / min						
Size Φ × L(mm):	12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
ISO No. : 658...	104 371 514 120	104 371 514 180	104 041 514 120	104 041 514 150	104 042 514 120	104 304 514 120	104 303 514 150



SK-2/4
Silicone
Polishers



Silicone Polishers
Special for resin base.
Material:silicon rubber mixing with normal abrasive.
Usage:For smoothing and initial polishing porcelain/
resin teeth.

Caution: Please use the products in allowed RPM.OR it
mayTake the products damaged!

Grit : Extra Coarse								
Packing(pcs): 12/100								
Order No. :								
Recommended R.P.M:	20,000-30,000 / min				15,000-20,000 / min			
Size Φ × L(mm):	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
ISO No. : 658...	104 243 544 040	104 243 544 060	104 300 544 100	104 275 544 100	104 237 544 100	104 292 544 130	104 012 544 140	104 145 544 150

Grit : Coarse								
Packing(pcs): 12/100								
Order No. :								
Recommended R.P.M:	20,000-30,000 / min				15,000-20,000 / min			
Size Φ × L(mm):	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
ISO No. : 658...	104 243 534 040	104 243 534 060	104 300 534 100	104 275 534 100	104 237 534 100	104 292 534 130	104 012 534 140	104 145 534 150

Grit : Mid-Coarse								
Packing(pcs): 12/100								
Order No. :								
Recommended R.P.M:	10000-20,000 / min							
Size Φ × L(mm):	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
ISO No. : 658...	104 243 524 040	104 243 524 060	104 300 524 100	104 275 524 100	104 237 524 100	104 292 524 130	104 012 524 140	104 145 524 150

Grit : Extra Fine								
Packing(pcs): 12/100								
Order No. :								
Recommended R.P.M:	10,000-15,000 / min				8,000-12,000 / min			
Size Φ × L(mm):	4 x 12	6 x 16	10 x 20	10 x 24	10 x 17	13 x 20	14 x 17	15 x 17
ISO No. : 658...	104 243 504 040	104 243 504 060	104 300 504 100	104 275 504 100	104 237 504 100	104 292 504 130	104 012 504 140	104 145 504 150

SG-0
Silicone
Polishers



Silicone Polishers
Special for precious metals wokpiece.
Material:silicon rubber mixing with normal abrasive.
Usage:For finishing and poishing precious/semi-precious
dental workpiece,such as finishing,smoothing and polishing
golden teeth and goleninner crown.

Caution: Please use the products in allowed RPM.OR it
mayTake the products damaged!

Grit : Coarse								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	8,000-10,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 658...	900 371 534 160	900 372 534 160	900 303 534 160	900 371 534 220	900 372 534 220	900 303 534 220	900 114 534 060	900 114 534 070

Grit : Mid-Coarse								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	8,000-10,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 658...	900 371 524 160	900 372 524 160	900 303 524 160	900 371 524 220	900 372 524 220	900 303 524 220	900 114 524 060	900 114 524 070

Grit : Fine								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	8,000-10,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 658...	900 371 514 160	900 372 514 160	900 303 514 160	900 371 514 220	900 372 514 220	900 303 514 220	900 114 514 060	900 114 514 070

Grit : Extra Fine								
Packing(pcs): 100								
Order No. :								
Recommended R.P.M:	8,000-10,000 / min							
Size Φ × L(mm):	16 x 1	16 x 2	16 x 2	22 x 1	22 x 3	22 x 3.5	6 x 23	7 x 20
ISO No. : 658...	900 371 504 160	900 372 504 160	900 303 504 160	900 371 504 220	900 372 504 220	900 303 504 220	900 114 504 060	900 114 504 070

Silicone
Polishers

Silicone Polishers
Special for precious metals wokpiece.
Material:silicon rubber mixing with normal abrasive.
Usage:For finishing and poishing precious dental workpiece,such as finishing,smoothing and polishing golden teeth and goleninner crown.
Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Silicone Polishers
Special for precious metals wokpiece.
Material:silicon rubber mixing with normal abrasive.
Usage:For finishing and poishing precious/semi-precious dental workpiece,such as finishing,smoothing and polishing golden teeth and goleninner crown.
Caution: Please use the products in allowed RPM.OR it mayTake the products damaged!

Silicone
Polishers

Grit : Coarse

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

							
SG 201 2	SG 202 2	SG 211 2	SG 212 2	SG 213 2	SG 216 2	SG 221 2	SG 222 2
10,000-12,000 / min							
3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
104 110 534 030	104 113 534 050	104 243 534 030	104 243 534 040	104 243 534 060	104 174 534 040	104 030 534 090	104 026 534 180

Grit : Coarse

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

						
SG 231 2	SG 235 2	SG 232 2	SG 234 2	SG 233 2	SG 241 2	SG 242 2
10,000-12,000 / min						
12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
104 371 534 120	104 371 534 180	104 041 534 120	104 041 534 150	104 042 534 120	104 304 534 120	104 303 534 150

Grit : Mid-Coarse

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

							
SG 201 3	SG 202 3	SG 211 3	SG 212 3	SG 213 3	SG 216 3	SG 221 3	SG 222 3
10,000-12,000 / min							
3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
104 110 524 030	104 113 524 050	104 243 524 030	104 243 524 040	104 243 524 060	104 174 524 040	104 030 524 090	104 026 524 180

Grit : Mid-Coarse

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

						
SG 231 3	SG 235 3	SG 232 3	SG 234 3	SG 233 3	SG 241 3	SG 242 3
10,000-12,000 / min						
12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
104 371 524 120	104 371 524 180	104 041 524 120	104 041 524 150	104 042 524 120	104 304 524 120	104 303 524 150

Grit : Fine

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

							
SG 201 4	SG 202 4	SG 211 4	SG 212 4	SG 213 4	SG 216 4	SG 221 4	SG 222 4
10,000-12,000 / min							
3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
104 110 514 030	104 113 514 050	104 243 514 030	104 243 514 040	104 243 514 060	104 174 514 040	104 030 514 090	104 026 514 180

Grit : Fine

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

						
SG 231 4	SG 235 4	SG 232 4	SG 234 4	SG 233 4	SG 241 4	SG 242 4
10,000-12,000 / min						
12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
104 371 514 120	104 371 514 180	104 041 514 120	104 041 514 150	104 042 514 120	104 304 514 120	104 303 514 150

Grit : Extra Fine

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

							
SG 201 5	SG 202 5	SG 211 5	SG 212 5	SG 213 5	SG 216 5	SG 221 5	SG 222 5
10,000-12,000 / min							
3 x 7	5 x 12	3 x 7	4 x 12	6 x 16	4 x 12	9 x 8.5	18 x 3
104 110 504 030	104 113 504 050	104 243 504 030	104 243 504 040	104 243 504 060	104 174 504 040	104 030 504 090	104 026 504 180

Grit : Extra Fine

Packing(pcs): 12/100

Order No. :

Recommended R.P.M:

Size $\Phi \times L(mm)$:

ISO No. : 658...

						
SG 231 5	SG 235 5	SG 232 5	SG 234 5	SG 233 5	SG 241 5	SG 242 5
10,000-12,000 / min						
12 x 1	18 x 1	12 x 2	15 x 2	12 x 4	12 x 2	15 x 2
104 371 504 120	104 371 504 180	104 041 504 120	104 041 504 150	104 042 504 120	104 304 504 120	104 303 504 150

Others



Product name :	Diamond Dresser
IDescription :	80*25*10mm
Order No :	DC001

EF0、EH0
Cutting Disc



粒度 Grit : 极细 Extra Fine			
包装 Packing(pcs): 6/12			
订货编号 Order No. :			
推荐转速 Recommended R.P.M:	12,000-20,000 / min		
规格 Size $\Phi \times L$ (mm):	38 x 0.6	22 x 0.2	38 x 0.8
ISO No. : 802...	900 371 514 380	900 371 514 220	900 371 514 380

