

CONTENTS

Dental Silicone Impression Material



Tangshan UMG Medical Instrument Co., Ltd.

ADD: No. 130-13 Torch Road, High-tech District, Tangshan City

TEL:+86-315-6314398 +86-186-33395498

FAX:+86-315-6314577

EMAIL:info@umgroupen.com

WEBSITE: https://www.umgroupen.com

Full way services and technical support offered by our professional engineers.

UMST-0



Description

heavy body

- 1. 0 contraction coefficient 2.no increase in occlusion
- 3. no rare delicacy
- 4. no injury to hands after long-term use
- 5.No need to wait, direct mold pouring
- 6.It takes only 70 seconds to solidify crown & Bridge.

Packing

2 spoons 262ml/barrel, 262ml*2/set, 12sets/ctn

UMST-2 dark blue



Description


medium body

- 1. Plastic stability, easy to operate
- 2. High precision, one-step in place

Packing

6pcs green mixing tips & 2 pcs white intraoral tips 50ml/pc,

2 pcs/box, 30 boxes/ctn



UMST-1 light blue




Description

heavy body

- 1. Plastic stability, easy to operate
- 2. High precision, one-step in place

Packing

6pcs green mixing tips & 2 pcs white intraoral tips 50ml/pc, 2 pcs/box, 30 boxes/ctn



UMST-2 green




Description

bite registration

- 1. High precision, low deformation rate
- 2. Fast curing in mouth and high hardness after curing
- 3. It is convenient to repair and has little bite resistance

Packing

6pcs green mixing tips & 2 pcs white intraoral tips 50ml/pc,



UMST-3 pink



Description


extra light body

- 1. High precision, good toughness, not easy to fracture
- 2. High hardness, not easy to deform
- 3. Good fluidity and fine detail impression

Packing

6pcs yellow mixing tips & 2 pcs yellow intraoral tips 50ml/pc,

2 pcs/box,30 boxes/ctn



UMST-3 green



Description

light body

- 1. High precision and good toughness
- 2. No deformation
- 3. Suitable for most mold taking

Packing

6pcs yellow mixing tips & 2 pcs yellow intraoral tips 50ml/pc,

4 pcs/box, 21 boxes/ctn



UMST-3 gray



Description

fit checker


- 1. Good thixotropy, easy to operate
- 2. The film thickness is small
- 3. Perfect permeability, convenient for observation

Packing

6pcs yellow mixing tips & 2 pcs yellow intraoral tips 50ml/pc, 2 pcs/box, 30 boxes/ctn

UM-YD33D printer

Everything can be printed







NO need to level ready to ues



Adjustable thickness Free conctrol



Quality is excellent on the safe side



Anneng temperature control Peace of mind printing



Double line guide rail Eliminate ripples



All-around printing Easy to do

Z-axis accuracy: 0.01mm

Printing accuracy: ±0.025mm

Printing speed: 75mm/h(layer thickness 0.1mm)

XY resolution: 2560× 1620mm

Print layer thickness: 0.025~0.2mm (adjustable)

Projection screen size: 152 ×96mm Light source: wavelength 405nm

Print area: 81× 128mm Print height: 150mm

Platform size: 82× 140mm

Connection: U disk, WiFi, high-speed Ethernet

Rated power: 110w

Rated frequency: 50/60hz

Equipment size: 250×290 × 420mm

Equipment weight: 12kg



Temporary Crown



Guide Plate



Dental Mold