

Iodine Volts Swab



Contains: The tube contains povidone iodine solution about 0.45%-0.5%

Product specification: length 8cm

Note: it is forbidden for allergic people to use this product. Do not let children use it alone



Features: mild wound disinfection low pain

Scope of application:

1. Disinfect the skin of the umbilical cord of newborns.
2. Disinfection of small skin wounds
3. Skin disinfection before injection and infusion



Iodophor disinfectant

+



Cotton swab

=



Providone Iodine Swabs

Traditional iodophor swab



After opening, it is easy to run iodine, and the longer the time, the weaker the disinfection effect.



Broken iodophor cotton stick




Broken iodophor individually packaged in sealed packaging, no need to worry about iodine escape.



Run iodine problem?

How to use the broken cotton swab



- Step 1: With the colored ring facing upwards, push the cotton swab out of the packaging film
- Step 2: Break off the position of the color ring and wait for the iodophor to flow to the lower end
- Step 3: Just wipe the disinfection position

Easy to carry around



Real shot

